

## SAFETY BASICS *Healthcare Series*

### Bathtub, Shower and Toileting

**Activity:** Resident Handling

**Description:** Shower and toileting chairs

**When to Use:** Showering and toileting residents who are partially dependent, have some weight bearing capacity, can sit up unaided, and are able to bend hips, knees and ankles.



**Points to Remember:** Ensure that wheels move easily and smoothly; chair is high enough to fit over toilet; chair has removable arms, adjustable footrests, safety belts, is heavy enough to be stable, and that the seat is comfortable, accommodates larger residents, and has a removable commode bucket for toileting.

Ensure that brakes lock and hold effectively on at least two wheels and that weight capacity is sufficient for resident handling.

Safety Basic Health Care Series Information is extracted directly from OSHA safety materials. For additional information consult [www.OSHA.gov](http://www.OSHA.gov)